RNase Remover

Code No. 9038

Size : 250 ml

Shipping at RT

Stored at RT

Lot No.

Contents

RNase Remover	250 ml	1 bottle
Spray bottle attachment		1 piece

Warning :

1. Wear plastic gloves when handling RNase Remover .
Contact of solution to skin can cause mild irritation and is harmful if swallowed. In case of contact with eyes, flush with water and contact physician. If swallowed, drink plenty of water and contact physician ; do not induce vomiting.
2. When using the spray bottle attachment , always work in a fume hood as inhalation of aerosol or vapor may cause irritation to lungs and mucous membranes.
3. Note that RNase Remover should not be used on corrodible metal surfaces.
4. Do not dilute because dilution will reduce its effectiveness.

Description : RNase Remover is a cleaning solution that is highly active against RNase contamination on the glass and plastic surfaces. It is possible to apply directly, and also to use with the spray bottle attachment.

Storage : Room temperature.

It is stable for 6 months at room temperature. It may generate a precipitate when stored at lower temperatures. The precipitate can be easily brought into solution by incubating at 37°C.

Instructions :

1. Decontamination of laboratory surfaces
Apply RNase Remover directly to the lab surface. Wipe thoroughly with a paper towel, rinse with sterilized distilled water and dry with a clean paper towel.
2. Decontamination of laboratory apparatus
Generously apply RNase Remover to a paper towel and wipe all exposed apparatus surfaces thoroughly. Rinse with a paper towel soaked with sterilized distilled water and dry with a clean paper towel. To clean small parts, briefly soak them in RNase Remover , rinse with sterilized distilled water well and dry with a paper towel.
3. Decontamination of plastic and glass vessels
Add RNase Remover enough to coat the entire surface of the vessel by swirling or vortexing. Discard the solution and rinse vessels thoroughly two times with sterilized distilled water.
4. Decontamination of pipettors
Following a manufacturer's instruction, remove the shaft from the pipettor and remove seals and gaskets from the shaft. Soak the shaft for one minute in RNase Remover , rinse the shaft thoroughly with sterilized distilled water, let dry and reassemble.

Note

For research use only. Not for use in diagnostic or therapeutic procedures.

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